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VITALITY HEALTH CHECK – VHC Vitamin D

Quantitative Vitamin D Test

Instructions for use


English

Version: 03/2022-06, 30.06.2022

VITALITY HEALTH CHECK - VHC Vitamin D

Quantitative vitamin D test

A Rapid "Sandwich" Immunochromatographic Test for the Quantitative Detection of total 25-OH Vitamin D in Human Whole Blood

For In Vitro Diagnostic use only

Read instructions for use before use

Intended Use

The VHC Vitamin-D Quantitative Vitamin D Test is an immunochromatography-based one step in vitro test. It is designed for the quantitative determination of total 25-hydroxy Vitamin D (25-OH Vitamin D) in human finger-prick blood. This assay provides a preliminary diagnostic test result and can be used for screening of Vitamin D deficiency. The liquid chromatography with tandem mass spectrometry (LC-MS/MS) assays or other quantitative immunoassays are recommended to further confirm the diagnostic test results.

Summary and explanation

Vitamin D is a steroid hormone responsible for enhancing intestinal absorption of calcium and the regulation of its homeostasis. The two common forms of Vitamin D are Vitamin D2 and Vitamin D3. Vitamin D3 is naturally produced in the human skin through the exposure to ultraviolet light and Vitamin D2 is mainly obtained from foods. Vitamin D is transported to the liver where it is metabolized to 25-hydroxy Vitamin D. In medicine, a 25-hydroxy Vitamin D blood test is used to determine Vitamin D concentration in the body. The blood concentration of 25-hydroxy Vitamin D is considered the best indicator of Vitamin D status. Vitamin D deficiency is now recognized as a global epidemic. Virtually every cell in our body has receptors for Vitamin D, meaning that they all require "Sufficient" Level of Vitamin D for adequate functioning. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin D deficiency has been linked to various serious diseases: Osteoporosis, Osteomalacia, Multiple Sclerosis, Cardiovascular Diseases, Pregnancy Complications, Diabetes, Depression, Strokes, Autoimmune Diseases, Flu, Different Cancers, Infectious Diseases, Alzheimer, Obesity and Higher Mortality etc. Therefore, now detecting (25- OH) Vitamin D level is considered as "Medically Necessary Screening Test", and maintaining sufficient levels not just to improve bone health, but to improve overall health and well-being. Multiple guidelines for Vitamin D deficiency have been published by various health organizations; but a common recommendation remained to be established. Recent literature (Worm 2010) has suggested a level classification for the classification of Vitamin D status (see table underneath). According to recent literature a level underneath 20 ng/ml is to be considered as critical low and it is recommended to contact a therapist for further diagnostics and treatment. Vitamin D levels between 40 and 60 ng/ml are regarded as optimal and have shown preventive effects. The multiple preventive effects of optimized Vitamin D are reflected by the health claims approved of the European Union regarding Vitamin D supplements: supporting conservation of normal Calcium levels, normal bones, normal muscle function, normal teeth, and normal function of the immune system.

Level	Vitamin D concentration range		Interpretation
1	< 20 ng/ml	< 50 nmol/l	Critically low vitamin D levels
2	21-30 ng/ml	52.5-75 nmol/l	Long-term vitamin D deficiency
3	31-40 ng/ml	77.5-100 nmol/l	Sufficient vitamin D level
4	41-60 ng/ml	102.5-150 nmol/l	Good vitamin D level
5	61-90 ng/ml	152.5-225 nmol/l	Very good vitamin D level
6	> 90 ng/ml	> 225 nmol/l	High vitamin D level

Note: Vitamin D concentrations are reported either in [ng/ml] or [nmol/l] (1 ng/ml is approximately 2.5 nmol/l). Make sure to check what unit system is used when comparing different test results!

VITALITY HEALTH CHECK – VHC Vitamin D

Quantitative Vitamin D Test

Instructions for use

English

Version: 03/2022-06, 30.06.2022

Test principle

VHC Vitamin-D Quantitative Vitamin D Test utilizes the principle of immunochromatography, a unique two-site “Sandwich” immunoassay on a membrane. The test employs a very “Exclusive” pair of anti-25-OH Vitamin D Monoclonal Antibodies; one conjugated with colloidal gold and another one immobilized on the solid phase. This will selectively detect Vitamin D with a high degree of sensitivity and specificity. As the test sample flows through the membrane assembly within the test device, the colored anti-25-OH Vitamin D-colloidal gold conjugate complexes with 25-OH Vitamin D from the sample. This complex moves further on the membrane by the capillary action to the test region (T) where it is immobilized by another anti-25-OH Vitamin D coated on the membrane, leading to formation of a pink / purple colored band, which confirms a positive test results. The intensity of colored band in the test line region is 25-OH Vitamin D concentration-dependent, higher the concentration of 25-OH Vitamin D in the tested sample, the stronger the colored band is. A control line is present in the test window to work as procedural control. This colored band should always appear on the control line region (C) if the test device is stored in good condition and the test is performed appropriately.

Materials provided

- Test cassettes for the quantitative vitamin D test VHC Vitamin-D (scope of the kit: 50 tests/pack, 25 tests/pack, 10 tests/pack)
- UniSampler sample collector (10 µl) and mixer including buffer (50 pieces in a pack of 50 tests, 25 pieces in a pack of 25 tests and 10 pieces in a pack of 10 tests)
- RFID card: 1
- Instruction for use: 1

Additional materials required:

- Timer or clock
- Safety lancet
- Alcohol swab
- VITALITY HEALTH CHECK Health Reader (VHC Reader) - (to be purchased separately)

Storage and stability

The test device should be stored at 4°C to 30°C and will be effective until the expiration date stated on the package. The product is humidity-sensitive and should be used immediately after being open. Any improperly sealed product should be discarded.

Precautions

1. For in vitro diagnostic use only.
2. Do not use the product beyond the expiry date.
3. Handle all specimens as potentially infectious.
4. The test cassette is a moisture-sensitive product, do not open foil pouch until it is ready to be tested.
5. Always compare the LOT number of the RFID card with the LOT number of the test batch, they must match.

Quality control


Good Laboratory Practice recommends the daily use of control materials to validate the reliability of device. If control values do not fall within established range, assay results are invalid. Control materials which are not provided with this test kit are commercially available. The VHC Vitamin-D Quantitative Vitamin D Test provides a built-in process control with a different antigen/antibody reaction at the control region (C). This control line should always appear regardless the presence of Vitamin D. If the control line does not appear, the test device should be discarded and the obtained result is invalid. The presence of this control band in the control region serve as 1) verification that sufficient volume is added, 2) that proper flow is obtained.

Sample collection and preparation


Note: The quantitative vitamin D test is designed for "decision point" blood (or serum) samples from the fingertip ONLY. **DO NOT use anticoagulated blood or plasma samples** for testing as anticoagulants may affect test results. The quantitative vitamin D test was developed for blood collection from the human finger. However, a vitamin D control or serum sample can also be used for the test. Instead of drawing blood from the fingertip with a blood collector, place 5µl of the vitamin D control sample or serum into the collection tube with a micropipette (not included in the kit) and follow the instructions for use.

- (1) Wash thoroughly by hand and dry completely.
- (2) Use the test at room temperature (21°C to 24°C)
- (3) Rub the finger selected for sampling slightly to increase blood flow.
- (4) Prick the side of the finger using a (safety) lancet.
- (5) Take 10 µl of blood with the capillary of the UniSampler sampler, mix with the sample buffer and perform the test immediately and read it after 15 minutes. The result after 15 minutes may not be accurate.


Instructions for using the UniSampler for sampling and sample application

- (1) The UniSampler consists of a Collection Tube filled with Sample Buffer and sealed (left) and a Blood collector with Cap (right).
- (2) Tap the Buffer tube on a Flat surface to ensure that all the liquid has come at the bottom.
- (3) Peel off aluminium foil seal from the top of the Collection Tube containing Sample Buffer.

- (4) Prepare the finger (clean and warm it) and prick the fingertip with the safety lancet.
- (5) Wait until a sufficiently large drop of blood has formed.
- (6) Gently touch the tip of Blood Collector to blood droplet. Capillary action will completely fill 10 µl blood and stop.
- (7) Check the blood level in the capillary of the Blood Collector. This must be completely filled (as shown above).

- (8) Fully insert Blood from the Blood collector into Collection Tube and push firmly to close tightly.
- (9) Shake the UniSampler with „Jerk“ 3 to 4 times to completely take out blood from Collector into the Sample Buffer, followed by complete mixing. **Mix the sample and the Buffer completely until a uniform colour is visible.**
- (10) Remove the cap of the UniSampler.
- (11) Invert the UniSampler and gently squeeze 3 drops of pre-mix blood into the Sample Well (S) of the Test cassette.

Important note on sample mixing and sample application – CAUTION!

- Complete (100%) PRE-MIXING of finger-prick blood with sample buffer is "EXTREMELY" important and CRITICAL Step to get correct result. This can be determined by checking the UNIFORM red color of pre-mix blood in Collection Tube and Blood Collector.
- Incomplete mixing of Blood with Buffer means Sample Preparation has been compromised, and the test result is likely to show lower values.
- Pressing of UniSampler should be "GENTLE" to get three full drops of pre-mix blood into the sample well (S).

Test implementation

- (1) Bring all materials and samples to room temperature.
- (2) Remove the test cassette from the sealed foil pouch and place it on a stable, flat surface.
- (3) Follow the instructions for using the UniSampler sampler and mixer.
- (4) After placing 3 drops of premixed blood in the sample well (S), read the result after 15 minutes using the VITALITY HEALTH CHECK Health Reader (VHC Reader).

Important note: The result may be inaccurate if the reading is taken after more than 15 minutes.


VITALITY HEALTH CHECK – VHC Vitamin D Quantitative Vitamin D Test

Instructions for use
Version: 03/2022-06, 30.06.2022

English

Quantitative determination with the VHC-READER

To perform the test, you need the test strip and the VHC-Reader (sold separately). For timing (15 minutes) you can use the VHC reader itself or your own timer.

Place the Adapter on top of the Test Device (test cassette) correctly.

Now place the VHC reader on the the top of the Adapter - make sure it is placed correctly.

Protocol without timer

Switch on the VHC reader by pressing the black button. The device now carries out a self-test, during which "WAIT" is displayed.

After an audible signal (Piep), "ON" is displayed. To take a measurement directly, press the button again ONCE FOR ONE SECOND. The display shows "RFID".

Place the LOT-specific RFID card included in the kit on the top of the VHC reader. This will upload the vitamin D test specific calibration data from the RFID card to the VHC reader. Following an audible signal (beep), "TEST" is displayed.


The meter now shows "RUN". After a successful data transfer, the measurement starts. The vitamin D concentration is displayed successively in ng/ml, in nmol/l and level 1 to 6.

Protocol with integrated timer

Switch on the VHC reader by pressing the black button. The unit now carries out a self-test, during which "WAIT" is displayed.

After an audible signal (Piep), "ON" is displayed. HOLD down the black button until "RFID" appears.

Place the LOT-specific RFID card included in the kit on the top of the VHC reader. This will upload the vitamin D test specific calibration data from the RFID card to the VHC reader. Following an audible signal (beep), "TEST" is displayed. Press the black button on the reader again and the 15-minute countdown of the timer will start.

The countdown display starts and after the 15 minutes have elapsed, the measurement is started automatically. The vitamin D concentration is displayed successively in ng/ml, in nmol/l and level 1 to 6.

VITALITY HEALTH CHECK – VHC Vitamin D


Quantitative Vitamin D Test

Instructions for use

English

Version: 03/2022-06, 30.06.2022

Standard curve of the VITALITY HEALTH CHECK READER (VHC Reader)

A typical standard curve for the VHC reader is shown on the left. The AU (area under curve) of the measured value is automatically converted into the units ng/ml and nmol/l in the VHC reader.

Performance features

Sensitivity

The sensitivity of the quantitative vitamin D test VHC Vitamin-D is 3 ng/ml (7.5 nmol/l). The sensitivity was determined by adding the mean value and 3.3 times the standard deviation for twenty vitamin D-free serum tests.

Detection range

The detection range of the quantitative vitamin D test VHC Vitamin-D with the VHC Reader is from 3 ng/ml (7.5 nmol/l) to 100 ng/ml (250 nmol/l).

Accuracy

The accuracy of the VHC vitamin D was determined by using human fingertip blood samples and comparing the results with those obtained from a reference ELISA with 25-OH vitamin D using the respective serum samples. The comparison showed a linear regression with a slope of 1.02 and a correlation coefficient of 92%. In conclusion, the results from the quantitative vitamin D test VHC vitamin D with human blood samples show good agreement with the results from the ELISA with the associated serum samples.

The accuracy of the quantitative vitamin D test VHC Vitamin-D was further evaluated using serum samples and compared to an LC-MS/MS assay ("gold standard" for the measurement of 25-OH vitamin D). The comparison showed a linear regression with a slope of 0.98 and a correlation coefficient of 98%. It can be deduced from this that the quantitative results from the VHC vitamin D strongly agree with the values obtained from the LC-MS/MS assay.

Precision

Sample	Number of Replicates	Mean (ng/ml)	Standard deviation	Coefficient of variation (CV)
Serum	20	40,4	3,1	7,6 %
Blood-1	10	33,0	2,5	7,6 %
Blood-2	10	38,7	3,9	10,1 %

Specificity

A total of 30 vitamin D-free serum samples were tested. All showed negative results, indicating a specificity of 100%.

No interferences and cross reactivities with Bilirubin, Triglycerides, Cholesterol, Vitamin B12 and Vitamin C were detected.

VITALITY HEALTH CHECK – VHC Vitamin D Quantitative Vitamin D Test

English

Instructions for use

Version: 03/2022-06, 30.06.2022




Expected results

VHC Vitamin-D Quantitative Vitamin D Test is a Rapid Quantitative assay. The test is intended to use for screening individuals to identify Vitamin D level. This assay provides only a preliminary analytical test result. The liquid chromatography with tandem mass spectrometry (LC-MS/MS) assays or quantitative immunoassays are recommended to confirm the analytical result.

Literature

- (1) Holick, MF. Vitamin D statuses: measurement, interpretation and clinical application. *Ann. Epidemiol.* 2009, 19(2):73-78.
- (2) Morris HA. Vitamin D: A Hormone for All Seasons - How much is enough? *Clin. Biochem. Rev.*, 2005, 26:21-32.
- (3) Moyad MA. Vitamin D: a rapid review. *Dermatol Nurs.* 2009, 21:25-30.
- (4) Zerwekh JE. Blood biomarkers of vitamin D status. *Am J. Clin Nutr.* 2008, 87:1087S-91S.
- (5) Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. *BMJ.* 2014, 348:g3656.
- (6) Worm N. Heikraft D: How the sun vitamin protects against heart attacks, cancer and other diseases of civilisation. *systemed Verlag, Lünen.* 2010, p12.

Symbols

	Article number
	Conformity marking
	Batch designation, LOT Number
	In vitro diagnostics
	Not for reuse
	Usable until (Expiry Date)
	Follow the instructions for use
	Content sufficient for x exams
	Manufacturer
	Temperature limit

Jungbrunnen - Fountain of Youth GmbH
Kleine Wasserstraße 23
D-18055 Rostock, Germany
info@jungbrunnenfoy.de
www.vitality-health-check.com

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test

English

Function and instructions for using the UniSampler for optimal sample preparation


Version: 03/2022-06, 30.06.2022

VITALITY HEALTH CHECK - VHC Vitamin D

Quantitative vitamin D test

An immunochromatographic "sandwich" rapid test for quantitative

Determination of total 25-OH vitamin D in human whole blood.

Function and instructions for using the UniSampler for optimal sample preparation - Training documents for professional users for self-study and reference

The VHC Vitamin D is an innovative rapid test that provides laboratory quality quantitative results within 15 minutes when used professionally. The procedure of the test within the cassette has been carefully developed to be robust and reliable. For the first time, the VHC Vitamin-D makes it possible to measure the important health parameter vitamin D level immediately and quantitatively, directly and everywhere at the point-of-care.

As a professional user in the healthcare sector, you are responsible for careful sample preparation and thus for the quality of the measurement. With the UniSampler you are able to be as precise as the machine. To support you in the use and processing of the VHC Vitamin D test, we have created this information. It serves as a training document for professional users.

It is very important that the user is concentrated during sampling and sample preparation and controls his steps. If the quantity of the sample to be applied is not correct or has not been mixed completely, this has a direct influence on the quality of the measurement result.

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test

English

Function and instructions for using the UniSampler for optimal sample preparation

Version: 03/2022-06, 30.06.2022

The UniSampler - your minilab for measuring, mixing and applying

For a quantitative test such as VHC vitamin D, sample preparation is crucial for the quality of the measurement. With the UniSampler you work at the point of care as precisely and controlled as in a professional laboratory. Therefore, it is important to know exactly which effects can occur that influence the sample preparation in the UniSampler.

The UniSampler is always actively involved in sampling and mixing:


- The capillary integrated in the sampler measures exactly the required 10 µl of capillary blood.
- For sample preparation, the capillary with the sample is placed on the sample container pre-filled with buffer so that the capillary is immersed directly.
- The sample is drawn from the capillary by **jerky movement and mixing** and mixed with the sample buffer.

The UniSampler was developed to optimally support you in measuring, mixing and applying. The integration of the sample capillary into the mixing vessel enables very precise work without aids if you use the UniSampler correctly and control the work steps.

The special role of the sample buffer for the test procedure

Vitamin D is fat-soluble. Therefore, it is also stored preferentially in fatty tissue. So that it can also be transported everywhere in the blood, it is bound there to the vitamin D-binding protein. The special composition of the sample buffer ensures that all vitamin D is released when the sample is mixed. This preparation step is important so that the specific antibodies reliably recognise and label all vitamin D in the sample.

- The sample buffer is applied. The vitamin D is free of vitamin D-binding protein.
- Continuous volume flow in the test strip. The antibodies in the test strip bind the vitamin D from the sample.
- The control band (C) forms first, the test line (T) rises dynamically until the end of the test (15 minutes).

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test

English


Function and instructions for using the UniSampler for optimal sample preparation

Version: 03/2022-06, 30.06.2022

Sampling with the UniSampler - measure 10 µl whole blood exactly!


It is important that the finger selected for sampling is warm and well supplied with blood in order to obtain a good sample quality. If necessary, it is recommended to immerse it in warm water for some time or to massage it lightly.

After puncturing with a suitable sterile lancet, often only a first small drop of blood comes out, which should be wiped away with an alcohol swab so that the wound does not close immediately. If the blood flow is only slight, the finger can be carefully massaged or lightly ring stitched but not "milked out". The second drop of blood that then emerges should be 2-3 mm in size to immediately fill the capillary completely. **Check this carefully!** If the capillary is not full or contains bubbles, this will affect the measurement. The mixing ratio of the sample is then not correct, which leads to incorrect measurements.

In case of doubt, it is always better to puncture a second time or use a second UniSampler capillary than to perform the VHC test with a bad sample!

When inserting the sampler into the sample collector of the UniSampler with the fully filled capillary pointing downwards until the tip of the capillary is immersed in the sample buffer, make sure that the capillary does not accidentally touch the rim or the inside of the sample collector.

Work carefully! Otherwise, part of the sample could be accidentally lost. Before mixing, close the UniSampler carefully by placing the cap with the sampler firmly on the sample collector until it fits tightly.

VITALITY HEALTH CHECK - VHC Vitamin D

English


Quantitative vitamin D test

Function and instructions for using the UniSampler for optimal sample preparation

Version: 03/2022-06, 30.06.2022

Sample and sample buffer in the UniSampler - Mix very carefully!

It is of particular importance for the quality of the mixture that sample and buffer are completely mixed. In this process, the sample in the capillary is ejected from the capillary by jerky movement with sufficient momentum and then completely mixed with the rest by swirling. Each user will optimise this mixing movement individually after a few tests. Work carefully and control the good sample mixing closely!

Only repeated, vigorous spinning and swivelling will result in complete mixing! Avoid formation of bubbles! Keep UniSampler tightly closed - Check sample for air bubbles in transmitted light.

In order to determine whether sample and sample buffer are completely mixed, you must look closely at the UniSampler several times in transmitted light. Initially, the capillary is still clearly more intensely coloured (A), because the capillary forces are holding back the sample. After complete emptying and mixing, the capillary finally appears somewhat lighter in transmitted light than sample mixture at the bottom of the UniSampler (C).

A

B

C

If you do not mix with sufficient momentum, a large part of the sample will remain in the capillary (A). If several drops form during careful mixing (B), unite them by slowly swirling the UniSampler. When completely mixed, the capillary appears lighter than the one at the bottom of the UniSampler (C).

Be sure to work carefully here and check the steps several times in transmitted light!

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test


English

Function and instructions for using the UniSampler for optimal sample preparation


Version: 03/2022-06, 30.06.2022

Sample application - Apply a full 3 drops to test cassette!

The lid of the UniSampler also contains the sample outlet or the drip tube with which the sample is applied. Here, too, it is important to measure correctly.

Bubble formation in the drip tube can lead to incorrect volumes. Make sure that no bubbles form here. Care must be taken when applying to the test cassette. If a drop touches a surface before it is fully formed, it will break off and then have too low a volume. Likewise, a trapped air bubble can change the volume.

**Work carefully! Drop three full drops on "S".
The volumes are important!**

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test

English

Function and instructions for using the UniSampler for optimal sample preparation


Version: 03/2022-06, 30.06.2022

Readout after exactly 15 minutes. The vitamin D signal increases dynamically!


After sample application, the continuous suction flow generated by the applied sample in the test cassette ensures controlled flow. Pay attention to the clear formation of the control line. The test line appears only after a few minutes and becomes continuously more intense.

So that the VHC reader "knows" exactly which measurement to take, each LOT (batch) of VHC vitamin D tests has its own RFID card with the same LOT numbers. Always make sure to check that you are using the correct RFID card. Compare the LOT numbers!

Position the adapter correctly on the test cassette! Otherwise the reader will take the measurement in the wrong place and the test result will not be correct!


Follow the steps of the test procedure.

Be sure to keep to the prescribed 15 minutes!

Check the test procedure, check the control band!

Test went well (C and T) Error images - Repeat test!

A B
A: strong signal;
B: Weaker signal;

C D E
C: test not passed (no T);
D: Control line (C) missing;
E: Control line (C) unclear;

VITALITY HEALTH CHECK - VHC Vitamin D Quantitative vitamin D test

English

Function and instructions for using the UniSampler for optimal sample preparation

Version: 03/2022-06, 30.06.2022

Some frequently asked user questions - FAQ

What should be done if the VHC reader indicates an error although the test has run correctly and the control band is clearly visible?

Check if you have accidentally inserted the test cassette the wrong way round in the adapter. In this case, the reader reads the control line or white area and reports an error or reads zero. Also check that the batteries in the VHC reader are still sufficiently charged. If necessary, you can also supply the VHC reader with power flexibly via the enclosed USB adapter cable.

What should be done if the control band is not clearly visible or streaky at the End of the test despite the sample being completely applied?

In this case, the test has not run properly inside the cassette and the test result must be discarded and repeated in any case. Please always take a photo of a test with a faulty control band. Please send this to us by e-mail (info@jungbrunnenfoy.de) together with the batch number (LOT) and your contact details and briefly describe what you have observed. We will get back to you as soon as possible.

What is the effect if the mixture of the sample buffer is not correct?

Without sufficient vigorous mixing, the sample will exit the capillary only slowly. If the emptying of the capillary is not precisely controlled, the mixture may contain too little sample. The measurement result will then deviate downwards.

What happens if less than 3 full drops of sample mixture are applied to the test strip?

The applied sample buffer volume creates a continuous suction flow in the strip of the test cassette, which drives the course of the test. However, if the sample inlet dries prematurely because the volume is too low, the test will not run fully for 15 minutes. The test line will then be too weak because the vitamin D signal will increase dynamically over the full 15 minutes. If you are unsure despite working carefully, apply another drop if necessary. It is important that the sample inlet does not run dry.

Why does the test have to be read out after 15 minutes?

The test is optimised for use at the point of care so that the result can be used quickly. If you follow the test procedure closely, you will see that the control band quickly becomes fully formed. The vitamin D test line, on the other hand, only appears with a delay and continues to rise dynamically the longer the test runs. It must be read out exactly after 15 minutes, because the calculation of the quantitative vitamin D value is based on the measured intensity of the test line.

How long do the batteries in the VHC reader last?

If you use the VHC reader exclusively for single measurements and not in timer mode, the batteries will last for up to 100 measurements. The battery requirement increases the more often you use the timer mode. If you mainly work in timer mode, you can also supply the VHC reader with power flexibly via the enclosed USB adapter cable.